

200 min Library Generation

- Mix 5 µl RNA, 5 µl **FS1x** and 5 µl **dT** **OR** 5 µl **Custom Targeted Primers**.
- Heat to 85 °C for 3 min, then cool to RT temperature. For targeted sequencing and FFPE samples skip this step.
- Pre-mix 4.5 µl **FS2x** and 0.5 µl **E1** per reaction. First Strand cDNA Synthesis
- Add 5 µl **FS2x** / **E1** mix per reaction.
- Incubate for 15 min at 37 °C (or up to 50 °C for targeted primers with a high Tm).
- Add 5 µl **RS1**, mix well. RNA Removal
- Incubate 10 min at 95 °C, then cool to 25 °C.
- Add 5 µl **RS2**, mix well.
- Pre-mix 13 µl **SS1x** and 2 µl **RSP** **OR** 2 µl **Custom Targeted Primers** per reaction, mix well.
- Add 15 µl **SS1x** / primer mix per reaction. Mix well. 2nd Strand Synthesis
- Incubate 1 min at 98 °C, slowly ramp down to anneal primers (25 °C - 50 °C).
- Incubate 30 min at 25 °C (or up to 50 °C for targeted primers with a high Tm).
- Pre-mix 4 µl **SS2** and 1 µl **E2** per reaction.
- Add 5 µl **SS2** / **E2** mix per reaction.
- Incubate 15 min at 25 °C (or 50 °C for targeted primers with a high Tm).
- Add 20 µl **PB** per reaction, mix well, incubate 5 min.
- Place on magnet for 2 - 5 min, remove and discard supernatant. Purification
- Add 40 µl **EB**, mix well, incubate 2 min at RT.
- Add 72 µl **PS** (or 48 µl **PS** for low input, FFPE, or low quality RNA), mix well, incubate 5 min at RT.
- Place on magnet for 2 - 5 min, remove and discard supernatant.
- Wash the beads (**PB**) twice with 120 µl 80 % EtOH, 30 sec.
- Dry beads for 5 - 10 min.
- Add 20 µl **EB**, mix well, incubate 2 min at RT.
- Place on magnet for 2 - 5 min, transfer 17 µl of the supernatant into a fresh PCR plate.

70 min Library Amplification

- Pre-mix 7 µl **PCR** and 1 µl **E3** per reaction, mix well.
- Add 8 µl **PCR** / **E3** premix to 17 µl of each purified library.
- Add 5 µl **BC** (from the 96-well plate) for each reaction, mix well. PCR
- PCR: 98 °C, 30 sec
 - 98 °C, 10 sec
 - 65 °C, 20 sec
 - 72 °C, 30 sec
 - 72 °C, 1 min; 10 °C, ∞

$$\left. \begin{array}{l} 98\text{ }^\circ\text{C, }10\text{ sec} \\ 65\text{ }^\circ\text{C, }20\text{ sec} \\ 72\text{ }^\circ\text{C, }30\text{ sec} \end{array} \right\} \times \text{cycle number}$$

use qPCR to determine cycle number

ATTENTION: Increase cycle number for low input RNA and samples with low target RNA content!
- Add 30 µl **PB** (or 27 µl **PB** for low quality or FFPE RNA) per reaction, mix well, incubate 5 min.
- Place on magnet for 2 - 5 min, remove and discard supernatant.
- Add 30 µl **EB**, mix well, incubate 2 min at RT.
- Add 30 µl **PS**, mix well, incubate 5 min at RT. Purification
- Place on magnet for 2 - 5 min, remove and discard supernatant.
- Wash the beads twice with 120 µl 80 % EtOH, 30 sec.
- Air dry beads for 5 - 10 minutes.
- Add 20 µl **EB**, mix well, incubate 5 min at RT.
- Place on magnet for 2 - 5 min, transfer 15 - 17 µl of the supernatant into a fresh PCR plate.

ATTENTION: Spin down solutions before opening tubes or plates!

The background of the page is white with a decorative pattern of semi-transparent blue spheres of various sizes. These spheres are connected by thin, light blue lines that create a network-like structure across the page. The spheres have a glossy, 3D effect with highlights and shadows.

QuantSeq-Flex Targeted RNA-Seq Library Prep Kit · Reference Card

Lexogen GmbH

Campus Vienna Biocenter 5

1030 Vienna, Austria

Telephone: +43 (0) 1 345 1212-41

Fax: +43 (0) 1 345 1212-99

E-mail: info@lexogen.com

© Lexogen GmbH, 2015